

# Characterization and evaluation of a novel C18 column based on organic/inorganic hybrid silica for HPLC and UHPLC

Takashi Sato, Noriko Shoji, Akiko Matsui, Mai Yamashita, Takatomo Takai and Naohiro Kuriyama  
YMC Co., Ltd., Ishikawa, Japan

## Introduction

Silica based reversed phase columns have been widely used for analytical and preparative chromatographic field. The silica based packing materials have low stability under alkaline conditions, and have a limited usable pH range. Recently, there has been much attention given to hybrid materials that have two aspects of inorganic and organic character to improve the chemical stability.

We have developed a new type of hybrid C18 stationary phase, YMC-Triart C18. Triart C18 is based on multi-layered organic/inorganic hybrid particles with 5 μm, 3 μm and a novel 1.9 μm diameter which are produced with a combination of our existing technologies for silica manufacturing and flow microreactor. We also have applied the optimized technology of surface modification for hybrid silica of Triart C18 to all particle sizes for improving durability, scalability, selectivity and peak shapes for various types of compounds.


In this poster, we will evaluate the pH stability and chromatographic performance of Triart C18 comparing commercially available columns, and show some application data utilizing characteristics of this material. And we would also present an example of method transfer from conventional HPLC using 5 μm particle to UHPLC using 1.9 μm particle.

## Features & benefits of YMC-Triart C18

- Three core technologies for particles and surface modification
  - 1. A multi-layered organic/inorganic hybrid particle
  - 2. A precise granulation with microreactor technology
  - 3. A proprietary C18 bonding and a multi-stage, multi-compound end-capping
- Symmetrical peak shapes and reproducible retention for all types of compounds under a variety of mobile phase conditions
- Improved speed and resolution in UHPLC analysis on 1.9 μm columns with operating pressure up to 100 MPa (14,500 psi)
- Superior column-to-column and lot-to-lot reproducibility provided by YMC's rigorous manufacturing control system
- Outstanding chemical and physical durability over a wide pH range at a high temperature


## Specification of YMC-Triart C18

Base material	Multi-layered organic/inorganic hybrid
Stationary phase	Polymerically bonded C18 group (USP class: L1)
Particle size	1.9 μm (New), 3 μm, 5 μm
Pore size	120 Å
Carbon loading	Approx. 20%
End-capping	Yes ("multi-stage end-capping" technology)
pH range	1-12
Temperature limit (Recommendation)	70°C for pH 1-7 50°C for pH 7-12



## Characterization and evaluation of YMC-Triart C18


### Comparison of column efficiency among different particle sizes



- 1.9 μm particle exhibits higher efficiency and maintains efficiency over a wide range of flow rate compared to 5 μm and 3 μm particles.
- 1.9 μm YMC-Triart C18 enables ultra high throughput analysis by using shorter length column and increasing flow rate.


### Comparison of selectivity and scalability among various ODS columns

#### Selectivity chart -Hydrophobicity and Hydrogen bonding capacity-




- YMC-Triart C18 has moderate hydrophobicity and hydrogen bonding capacity as a standard C18 phase and its selectivity is consistent across all of the particle sizes.
- Many of sub-2 and sub-3 μm UHPLC columns (including fully porous or superficially porous; shown in green) tend to have lower hydrophobicity than the common conventional HPLC columns (shown in black), and the selectivity of some phases varies between different particle sizes even in the same brand. It would limit a seamless method transfer between HPLC and UHPLC.

### Comparison of chromatographic performance for ionic compounds with UHPLC columns




### Evaluation of chemical durability under pH 9.5 condition

#### Analysis of barbiturates



- No change in retention time of barbiturates was observed even after 1000 runs at elevated pH.
- High chemical durability of YMC-Triart C18 achieved by applying hybrid particles and novel surface modification allows to utilize a wide pH range for better method development.

### Comparison of scalability across different particle sizes




- YMC-Triart C18 columns show the identical selectivity and the excellent peak shapes of basic compounds across all of the particle sizes including 1.9 μm. It allows predictable scale up from UHPLC to conventional HPLC and even to semi-preparative LC, and vice versa.
- There are some differences in selectivity, retention, and also peak shapes across the different particle sizes of commercially available C18 phases in the same brand (or alternative recommended by its manufacture).
- The Core-shell type C18 columns show significant peak tailing and have limited scalability because of lack of larger particle sizes.

## Specification of YMC-Triart C18

## Application


### Development of Ultra Fast separation method for 9 pesticides

#### Structures of pesticides




- The pesticides shown in left are used commonly on golf course.
- A large difference in chemical properties of those pesticides results in long analysis time for a precise separation.


### Increase throughput for the analysis of pesticides



### Lot-to-lot reproducibility of 1.9 μm YMC-Triart C18



### Comparison of peak shapes and selectivity among commercial UHPLC columns



## Conclusions

- The enhanced durability and chromatographic performance of YMC-Triart C18 offers the maximum flexibility in separation conditions across an expanded pH range.
- 1.9 μm YMC-Triart C18 with excellent chromatographic performance and 100 MPa of maximum operating pressure enables ultra-fast and reliable analysis.
- Identical chromatographic performance and selectivity of Triart C18 across different particle sizes provides mutual method transfer among UHPLC, HPLC and even semi-preparative LC.