

# DAD-50

**DAD** series

For medicine, chemistry, and food preparation

## **Preparative Dynamic Axial Compression Column**

Model DAD-50 is an automatic packing column developed for high-pressure preparation. Constant pressurization of the packed bed enables self-packing of a high-durability column without creating voids. The slurry port enables operation by the closed system from supplying of gel to pressurized packing with ease.

Higher sealing performance is achieved by a uniquely developed sealing structure (Japanese Patent No. 3929407). It can be used in a packing cycle for long periods. This column is entirely free from electric parts and therefore, safety use in explosion-proof area is possible. The column, manufactured in Japan, corresponds to IQ/OQ validation. We ensure a thorough support system for the maintenance and supply of consumables.

It is possible to select from various packing materials depending on the preparative object. Preparation and purification will be realized to your satisfaction together with a group of products designed exclusively by YMC.

Main Characteristics		
Product number	DAD-50-850S	
Inner diameter	50 mm	
Packed bed length	100 - 500 mm	
Maximum pressure	10 MPa	
Compression	Constant pressurization (DAC)	
Driving direction	Downward *	
Temperature	10 - 40 deg C, Below 50% RH	
Installation height	Over 3.0 m (Hoist is required)	
Dimensions	720(W) × 620(D) × 2560(H) mm	
Weight	About 80 kg	
Utility	0.5 MPa dry air	
Maximum Packing material	500 g	

Recommended packing material:

Silica, Organic/inorganic hybrid silica 10 µm and larger

Resistance to solvent:

Ordinary solvents used in reverse phase and normal phase HPLC

\* It is easy to unpack and collect gels from the bottom of column by moving the cylinder down.



### **Specifications**

Nai	me	Preparative Dynamic Axial Compression Column DAD-50
Product	number	DAD-50-850S
	Size	Inner diameter = $\phi$ 50 mm, Length = 850 mm
Body	Inner finishing	Electrolytic polishing after honing process
	Roughness	Ra ≤ 0.5
Filter	Material	SUS316L (Sintered filter)
Filter	Frit size	2 μm
	Material	PTFE
Seal	Туре	Uniquely developed sealing structure
Seal	Feature	Packing seals effectively in response to inner pressure of column.
		Easy replacement
	Source	Air-driven hydraulic pump
Cylinder	Туре	All stainless cylinder
	Required air	MAX 700 NL/min
	Switching valve	Change orientation of cylinder alignment (UP, DOWN, HOLD)
Operation	Pressure regulator	Adjust pressure of oil pump
	Pressure meter	Monitoring air operation and oil pressure
Wetted	material	SUS316, SUS316L, PTFE, PEEK
	Inlet/Outlet	1/8" Swagelok tube fitting
Connection	Slurry port	ISO 15A sanitary ferrule
	Air	$\phi$ 6 one-touch fitting
Funnel	Material/Fitting	SUS316 / ISO15A sanitary ferrule

#### **Option**

Slurry container	Feature	Slurry can be carried in closed system.
	Product number	SC-2-20
	Pump	Air-driven diaphragm pump
	Required air	MAX 200 NL/min
	Wetted material	SUS316, SUS316L, PTFE, PFA, FEP
	Inner finishing	Electrolytic polishing
	Connection	Column connection: ISO 15A sanitary ferrule
		Air connection: $\phi$ 6 one-touch fittings
Cylind	der oil	USDA, NSF H1 approved oil is available.

For more information or to place an order, please contact us. Please feel free to ask our specialists about products or service.



www.ymc-group.com



sales@ymc.co.jp

#### **Worldwide Availability**

YMC America, Inc. www.ymcamerica.com

YMC Europe GmbH www.ymc.de

YMC India Pvt. Ltd. www.ymcindia.com

YMC Korea Co., Ltd. www.ymckorea.com

YMC Shanghai Rep. Office www.ymcchina.com

YMC Singapore Tradelinks Pte. Ltd. www.ymc.sg

YMC Taiwan Co., Ltd. www.ymctaiwan.com



#### YMC CO., LTD.

YMC Karasuma-Gojo Bldg., 284 Daigo-cho Karasuma Nishiiru Gojo-dori, Shimogyo-ku Kyoto, 600-8106, Japan

TEL: +81-75-342-4515 FAX: +81-75-342-4550 www.ymc.co.jp sales@ymc.co.jp

The above information is as of December, 2014 and is subject to change without notice. ©2014 YMC CO., LTD. All Rights Reserved.

Distributor