

Care and Use Instructions Ion Exchange Media BioPro IEX SmartSep Q / BioPro IEX SmartSep S

1. Introduction

Thank you for purchasing YMC products for ion exchange media. BioPro IEX SmartSep Q/S are based on high mechanical strength hydrophilic polymer beads with a strong-anion exchanger (quaternary ammonium group) or a strong-cation exchanger (sulfo group), and are most suitable for biopharmaceutical purification.

BioPro IEX SmartSep Q/S, which are manufactured under highly controlled conditions, must pass a series of stringent tests before being accepted for shipment (Please refer to the certificate of analysis). To ensure optimal performance and durability of the media, please follow these instructions.

2. Specification

| Item | Strong-anion exchanger BioPro IEX SmartSep Q | | | Strong-cation exchanger BioPro IEX SmartSep S | | |
|------------------|--|---|----|---|---|----|
| | Particle size (µm) | 10 | 20 | 30 | 10 | 20 |
| Matrix | Hydrophilic polymer beads | | | | | |
| Functional group | -R-N ⁺ (CH ₃) ₃ | | | -R-SO ₃ ⁻ | | |
| pH range | 2.0 – 12.0 | | | 2.0 – 12.0 | | |
| Temp. range (°C) | 4 – 40 | | | 4 – 40 | | |
| Pressure range | Regular use: 3 MPa or less than 3 MPa Upper limit: 4 MPa | Regular use: 2 MPa or less than 2 MPa Upper limit: 3 MPa | | Regular use: 3 MPa or less than 3 MPa Upper limit: 4 MPa | Regular use: 2 MPa or less than 2 MPa Upper limit: 3 MPa | |
| Shipping solvent | 20% ethanol aqueous solution | | | | | |
| Application | Highly efficient polishing step and intermediate purification step | | | | | |

3. Packing instructions

3-1 Resin preparation (Removal of fines)

Recommended slurry solvent : Deionized water, buffer solution or 20% ethanol solution

Amount of media required : approximately 1.05 X to 1.10 X column volume

1. Suspend the settled resin in the shipping container(s) and transfer the slurry into a clean container with a capacity of 5 times or larger than the amount of resin.
2. Allow the slurry to settle for a certain time. Then measure the amount of the resin and decant the supernatant.
3. Add four times the resin volume of the slurry solvent stated above to the container.
4. Stir the solution with a rod or paddle to make homogeneous slurry. Do not use a sharp-edged paddle or a magnetic stirrer to avoid fine generation.
5. Allow the 30 µm resin to settle for 120 minutes. The smaller the particle size is, the more time for sedimentation is required.
6. Once resin has settled, decant the supernatant.
7. Repeat Step 3-6 until the supernatant at the Step 6 becomes clear.

3-2 Packing slurry preparation and column packing

Recommended packing solvent : A solution with high ionic strength (a component of mobile phase which has the highest ionic strength. For example, 1 M NaCl or 0.5 M Na₂SO₄)

1. Filter the slurry solution obtained after 3-1. If 20% ethanol is used as the slurry solvent at 3-1, rinse the resin by filtration with deionized water (This process can be skipped if deionized water or buffer is used at 3-1).
2. Rinse the resin on the filter with approximately three times the resin volume of the packing solvent stated above.
3. Transfer the resin to a clean container with a capacity of 5 times or larger than the amount of resin.
4. Add the packing solvent to adjust the slurry concentration to 30-50%. The slurry concentration can be calculated with a following equation. Slurry concentration (%) = Resin volume (L) / Total slurry volume (L) X 100
5. Stir the solution with a rod or paddle to make homogeneous slurry.
6. Once homogeneous slurry is formed, gently pour it into the column. Prevent air coming into the column during transfer.
7. Connect a pump to the column and start pumping the packing solvent at approximately two times of the actual flow rate in use. Maintain the flow rate until the packed bed becomes stable. Then disconnect the pump. The desired column volume can be obtained by flowing the packing solvent through the column and compressing the resin.
8. When using an axial compression column, the column performance may be improved by pushing in after pumping..

NOTE) Also refer to the Instruction manual of the column using and follow the packing procedures described in it, when necessary.

3-3 Column performance evaluation

When the column packing is completed, evaluate the column performance evaluation by injecting a sample to the column, and determine the column theoretical plate number (N/m) and peak asymmetry factor (As). Typical evaluation condition is shown below. If there is a large difference between the value(s) obtained and typical performance listed below, adjust the packing parameter(s) (slurry concentration, flow rate of the slurry solvent etc.) and repeat the packing procedure.

Example of test conditions for the column packing evaluation

| Condition 1 | | Condition 2 | | Column performance | | | |
|------------------|--------------------------------|--|--|--------------------------------|---------------|---------------|--------------|
| Mobile phase | 0.5 M NaCl | Low ionic strength buffer | | Particle size | 10 μ m | 20 μ m | 30 μ m |
| | | Strong anion exchanger (Q): 20 mM Tris-HCl buffer (pH 8) | | Theoretical plate number (N/m) | \geq 20,000 | \geq 10,000 | \geq 7,000 |
| | | Strong cation exchanger (S): 20 mM Phosphate Buffer (pH 7) | | | | | |
| Sample | 1 M NaCl | Formamide (2 μ L/mL) | | Asymmetry factor (As) | 0.7~1.4 | | |
| Detection | Conductivity | UV at 220 nm | | | | | |
| Flow rate | Approximately 70 to 90 cm/hour | | | | | | |
| Temperature | Room temperature (25°C) | | | | | | |
| Injection volume | 1% of bed volume | | | | | | |

* These are just some guide values. Desired separation could be achieved even if the values obtained are out of the range shown above, depending on your usage and/or application.

* Sample diffusion in the system flow path (extra-column volume) greatly affects column performance. If the packing conditions do not change the column performance, check if the inner diameter of connecting tubes and equipment performance are appropriate.

4. Equilibration and elution

- Equilibrate with about 5–10 column volumes of initial mobile phase before using a column for chromatographic separations.
- Generally, samples are adsorbed onto the top of the column with 20–50 mM of buffer as the initial mobile phase, and then eluted with a salt-concentration gradient method (for example, adjusting sodium chloride concentration between 0 and 0.5 M) or pH gradient method. It is recommended to flush the column with buffer containing about 1 M of sodium chloride after each run in order to remove residual impurities in the column staying with the final mobile phase.
- Water-soluble organic solvent (maximum of 30%), can be added to the mobile phase. Before adding such solvent, make sure that salt in the buffer will not precipitate.

5. Cleaning and regeneration

A change of retention time or peak shape and/or pressure increase might result from the adsorption of fat-soluble substances or precipitated impurities in a sample. In such cases, follow these steps for column cleaning and regeneration. If these procedures will not solve the problem, we recommend using new media.

5-1 Common cleaning methods

- Cleaning in place (CIP) - The following CIP is effective when there is a change in column performance or before long-term storage (it is recommended to wash the column without connecting to a detector).

First, 3 to 5 CV (CV : Column Volume) of 1 to 2 M NaCl is flowed through the column. Then, 3 to 5 CV of 0.1 to 0.5 M NaOH is flowed through the column. The cleaning efficiency of the column can be improved by increasing the concentration of NaOH (up to 1 M) and/or the exposure time by slowing the flow rate. In order to neutralize the column, 3 to 5 CV of 1 to 2 M NaCl is flowed through the column. After neutralization, equilibrate the column thoroughly with the same mobile phase used for the next step.

The state of contamination or type of washing solution can cause a pressure increase. In such cases, reduce the flow rate.

- Batch method

Soak and agitate the media in washing solution about 3–5 times of the media volume. After leaving it stationary, remove the supernatant fluid by decantation. Repeat the process 2–3 times. As a washing solution, the solvents used for CIP can be used.

5-2 Cleaning with surfactants and other additives

- Other additives such as urea (\leq 8 M) or guanidine hydrochloride (\leq 6 M), which are commonly used as protein denaturants are useful. Avoid solvents containing oxidant for the mobile phase.
- Nonionic surfactants, cationic surfactants or anionic surfactants are useful, but avoid anionic surfactants for BioPro IEX SmartSep Q and avoid cationic surfactants for BioPro IEX SmartSep S.
- The state of contamination or type of washing solution (high viscosity, etc.) can cause a pressure increase. In such cases, reduce the flow rate.

6. Storage

Store the bulk media in 20% ethanol aqueous solution. Store the bulk media in the original container at temperatures from 4 to 35°C. Keep the container closed tightly.