



# YMC Oligonucleotide Columns




IP-RP  
AEX  
LC/MS





## Content/Phase selection guide

|                                         | Page  |
|-----------------------------------------|-------|
| Oligonucleotide Application Notes ..... | 4–8   |
| IP-RP Oligonucleotide Columns .....     | 10–15 |
| AEX Oligonucleotide Columns .....       | 18–23 |
| Ordering Information.....               | 26–30 |





### Oligonucleotides / Nucleic Acids






### Oligonucleotide Applications




## Oligonucleotide applications

## Oligonucleotides d(T)2-20 method transfer from HPLC to UHPLC

Conventional LC method




UHPLC method



Eluent: A) 10 mM di-n-butylamine-acetic acid (pH 6.0)  
B) methanol  
Temperature: 37°C

Detection: UV at 269 nm  
Injection: 1  $\mu$ L (5 nmol/mL)

## Synthetic oligonucleotides




Column: YMC-Triart C18 (1.9  $\mu$ m, 12 nm) 50 x 2.1 mm ID  
Part No.: TA12SP9-05Q1PT  
Eluent: A) 200 mM HFIP\*–8 mM triethylamine  
B) methanol  
Gradient: 10–20%B (0–10 min)

Flow rate: 0.42 mL/min  
Temperature: 65°C  
Detection: UV at 260 nm  
Injection: 1  $\mu$ L (2–4 nmol/mL)  
Courtesy of M. Yamada, SHIMADZU CORPORATION (JP)

\*hexafluoroisopropanol

## LC-MS analysis of synthetic 27-30 mer oligonucleotides

Sample: Primer of DNA sequencing  
5'-CCGCTCGAGCTAAAAAAGCCTGTGTACC-3' (30 mer)




Column: Hydrosphere C18 (3  $\mu$ m) 50 x 2.0 mm ID  
Part No.: HS12S03-0502WT  
Eluent: A) 10 mM DBAA (pH 6.0)  
B) Mobile phase A / acetonitrile (50/50)  
Gradient: 58%–62%B (0–20 min), 62%B (20–25 min)

Flow rate: 0.2 mL/min  
Temperature: 35°C  
Detection: UV at 269 nm and ESI negative-mode  
Injection: 1  $\mu$ L (10 pmol/component)



## LC/MS analysis of miRNA

RP



Courtesy of M. Yamada, SHIMADZU CORPORATION

5'-pUGG AGU GUG ACA AUG GUG UUG-3'  
5'-pUGG AGU GUG ACA AUG GUG UUG U-3'


(21 nt, MW 6,890.1)  
(22 nt, MW 7,196.3)

Column: YMC-Triart C18 (3  $\mu$ m, 12 nm) 150 x 2.0 mm ID  
Part No.: TA12S03-1502WT  
Eluent: A) 10 mM di-n-butylamine-acetic acid (pH 7.5)  
B) 10 mM di-n-butylamine-acetic acid (pH 7.5)/acetonitrile (50/50)  
Gradient: 62–72%B (0–20 min)  
Flow rate: 0.2 mL/min

Temperature: 30°C  
Detection: A) UV at 260 nm  
B) ESI-negative mode  
Injection: 4  $\mu$ L (5 nmol/mL)  
System: LC) Shimadzu Prominence  
MS) Shimadzu LCMS2020

## Influences of mobile phase conditions on intensity of ESI-MS

RP




Column: Hydrosphere C18 (3  $\mu$ m) 50 x 2.0 mm  
Part No.: HS12S03-0502WT  
Flow rate: 0.2 mL/min  
Temperature: 35°C  
Detection: ESI negative mode  
Injection: 5  $\mu$ L (25 pmol/component)  
a) Eluent: A) 10 mM TEAA buffer (pH 6.0)  
B) eluent A/acetonitrile (80/20)  
Gradient: 50–65%B (0–20 min)  
b/c) Eluent: A) 10/5 mM DBAA buffer (pH 6.0)  
B) eluent A/acetonitrile (50/50)  
Gradient: 30–75%B (0–20 min)

d) Eluent: A) 10 mM TEAA buffer (pH 6.0)  
B) eluent A/methanol (50/50)  
Gradient: 44–50%B (0–25 min)  
e/f) Eluent: A) 10/5 mM DBAA buffer (pH 6.0)  
B) eluent A/methanol (20/80)  
Gradient: 42–70%B (0–20 min)  
g) Eluent: A) 400 mM HFIP-TEA buffer (pH 7.0)  
B) methanol  
Gradient: 7–35%B (0–20 min)  
Sample: Oligodeoxythymidylic acid [d(pT)<sub>2-20</sub>]

## Oligonucleotide applications

## LC-HRMS analysis of the antisense oligonucleotide Mipomersen (Kynamro®)

RP



Column: YMC-Triart C8 **metal-free** (1.9  $\mu$ m, 12 nm)<sup>1</sup> 100 x 2.1 mm ID  
 Part No.: TO12SP9-10Q1PTP  
 Eluent: A) water/triethylamine/HFIP<sup>2</sup>  
 (100/0.4/2; triethylamine 28.0 mM, HFIP 135.8 mM)  
 B) methanol/triethylamine/HFIP (100/0.4/2)  
 Gradient: [Sample separation step]  
 10–40%B (0–5.0 min)


[Column wash steps]  
 40–70%B (5.0–5.1 min), 70%B (5.1–7.0 min), 70–10%B (7.0–7.1 min),  
 10%B (7.1–8.0 min), 10–90%B (8.0–8.1 min), 90%B (8.1–9.0 min),  
 90–10%B (9.0–9.1 min), 10%B (9.1–10.0 min),  
 10–90%B (10.0–10.1 min), 90%B (10.1–11.0 min),  
 90–10%B (11.0–11.1 min)  
 Flow rate: 0.3 mL/min  
 Temperature: 50°C  
 Injection: 10  $\mu$ L (1000 ng/mL)  
 System: LC) Vanquish Binary Pump H system  
 HRMS) Orbitrap HRMS Q Exactive Plus

\*1 Prewash the column prior to the first use with water/methanol/phosphoric acid (70/30/0.1) for 1 hour  
 \*2 1,1,1,3,3,3-hexafluoro-2-propanol

Reference: Y. Sun et al, Development of a bioanalytical method for an antisense therapeutic using high-resolution mass spectrometry, Bioanalysis, 2020 NOV 26, doi: 10.4155/bio-2020-0225.

## Challenging phosphorothioate oligonucleotides

RP




Column: YMC-Triart C8 metal-free (1.9  $\mu$ m, 12 nm) 100 x 2.1 mm ID  
 Part No.: TO12SP9-10Q1PTP  
 Eluent: A) 15 mM triethylamine-400 mM HFIP<sup>△</sup>  
 B) methanol  
 Gradient: 10–20%B (0–20 min)

Flow rate: 0.42 mL/min  
 Temperature: 70°C  
 Detection: UV at 260 nm  
 Injection: 1  $\mu$ L (each 1.25 nmol/mL)

\*1,1,1,3,3,3-hexafluoro-2-propanol

## Separation of oligonucleotides modified with disulfides

RP



Column: 5  $\mu$ m, 250 x 4.6 mm ID  
Product code: TB30S05-2546PTH  
HS12S05-2546WT  
Eluent: A) 50 mM TEAA\* (pH 7.0)/acetonitrile (95/5)  
B) acetonitrile  
Gradient: 5–95%B (0–30 min), 95%B (30–35 min),  
95–5%B (35–35.1 min), 5%B (35.1–45 min)

Flow rate: 1 mL/min  
Temperature: 50°C  
Detection: UV at 260 nm  
Sample: Crude reaction mixture  
By courtesy of Saki Kawaguchi,  
Chemistry Department, Nagoya University, Japan

\*Triethylammonium acetate


## Easy purification of oligonucleotides with YMC-Actus semi prep columns

RP

### Purification of synthetic 30mer oligonucleotide

**Analysis** 1.0 mL/min, 5  $\mu$ L injection  
**Hydrosphere C18**  
50 x 4.6 mm ID, 5  $\mu$ m

**Purification** 19 mL/min, 100  $\mu$ L injection  
**YMC-Actus Hydrosphere C18**  
50 x 20 mm ID, 5  $\mu$ m




Part Nos.: HS12S05-0546WT  
HS12S05-0520WX  
Eluent: A) 10 mM DBA-acetic acid (pH 6.0) / methanol (60/40)  
B) 10 mM DBA-acetic acid (pH 6.0) / methanol (20/80)  
Gradient: 10%–35%B (0–30 min.)  
Temperature: ambient  
Detection: UV at 269 nm  
Sample: synthetic oligonucleotide (100  $\mu$ M)

**purity>99%**


# Oligonucleotide applications

AEX


## Separation of synthetic oligonucleotides (single-strand DNA)



## Separation of synthetic oligonucleotides (single-strand RNA)



## Separation of synthetic 2'OMe oligonucleotides (single-strand RNA)



Column: BioPro IEX QF (5  $\mu$ m) 100 x 4.6 mm ID  
Part No.: QF00S05-1046WP  
Eluent: A) 10 mM NaOH  
B) 10 mM NaOH containing 1.0 M NaClO<sub>4</sub>  
Gradient: 25–55%B (0–15 min), 100%B (15–20 min)

Flow rate: 1.0 mL/min  
Temperature: 25°C  
Detection: UV at 260 nm  
Injection: 4  $\mu$ L (5 nmol/L)

\*1,1,1,3,3,3-hexafluoro-2-propanol





IP-RP



## IP-RP – UHPLC/HPLC Selectivities


### Features

- Selection of C18, C8 and C4 columns
- For ion-pairing (IP-RP) UHPLC and HPLC
- pH- and temperature stable phases
- Superior reproducibility
- Optional bioinert hardware

### Selectivities for oligonucleotides


|                                                 | YMC-Triart C18               | YMC-Triart Bio C18           | Hydrosphere C18 | YMC-Triart C8                | YMC-Triart Bio C4            |
|-------------------------------------------------|------------------------------|------------------------------|-----------------|------------------------------|------------------------------|
| <b>Modification</b>                             | C18 (USP L1)                 | C18 (USP L1)                 | C18 (USP L1)    | C8 (USP L7)                  | C4 (USP L26)                 |
| <b>Particle Size / <math>\mu\text{m}</math></b> | 1.9, 3, 5                    | 1.9, 3, 5                    | 2, 3, 5         | 1.9, 3, 5                    | 1.9, 3, 5                    |
| <b>Pore Size / nm</b>                           | 12                           | 30                           | 12              | 12                           | 30                           |
| <b>pH range</b>                                 | 1.0–12.0                     | 1.0–12.0                     | 2.0–8.0         | 1.0–12.0                     | 1.0–10.0                     |
| <b>Temperature range</b>                        | pH < 7: 90°C<br>pH > 7: 50°C | pH < 7: 90°C<br>pH > 7: 50°C | 50°C            | pH < 7: 90°C<br>pH > 7: 50°C | pH < 7: 90°C<br>pH > 7: 50°C |

*Biocompatible  
hardware available!*



YMC-Triart **metal-free** columns are available for improved sensitivity and peak shape of coordinating compounds such as nucleotides or oligonucleotides.

## Enhanced retention and resolution even at low buffer concentrations




The separation of oligo(deoxythymidylic acids), d(pT)2-20, was compared using 100 mM or 10 mM triethylammonium acetate (TEAA) buffer, under the same gradient conditions. Both Hydrosphere C18 and YMC-Triart C18 showed enhanced retention and resolution compared to other commercially available C18 phases designed for oligonucleotide analysis, even at the low ion-pairing buffer concentration such as 10 mM.



The higher concentration provides stronger retention and superior resolution of oligonucleotides, although a lower concentration has the advantages of increasing the signal intensity and reducing system contamination in HPLC-MS analysis.

## IP-RP – Temperature and mobile phase effects


## Effect of mobile phase and column temperature on separation of siRNA duplex



## Versatile wide pH stability



## Metal-free column hardware ideal for oligonucleotide analysis




### Specifications

|                   |                                                                                    |
|-------------------|------------------------------------------------------------------------------------|
| YMC-Triart Phases | C18, Bio C18, C8, Bio C4                                                           |
| Particle Size     | 1.9, 3, 5 $\mu\text{m}$                                                            |
| Inner layer       | PEEK                                                                               |
| Outer layer       | Stainless steel                                                                    |
| Frit              | PEEK                                                                               |
| Pressure limit    | 1.9 $\mu\text{m}$ : 100 MPa (15,000 psi)<br>3/5 $\mu\text{m}$ : 45 MPa (6,525 psi) |


Special column connectors required.  
See ordering information recommendations.

## Improved sensitivity for coordination compounds




1. ATP



2. ADP



3. AMP




Column: YMC-Triart C18 (3  $\mu\text{m}$ ) 50 x 2.1 mm ID  
Part Nos.: TA12S03-05Q1PTP (**metal-free**) or TA12S03-05Q1PTH (regular hardware)  
Eluent: 5 mM  $\text{HCOONH}_4$   
Flow rate: 0.21 mL/min  
Temperature: 25°C  
Detection: UV at 265 nm  
Injection: 1  $\mu\text{L}$  (10 mg/mL)  
System: bioinert/"non-metal" HPLC system

Metal coordinating compounds, which have a phosphate group in their structure, tend to show poor peak shape due to interactions with metals, such as the stainless steel in column bodies and frits. By using the **metal-free** column hardware, better peak shapes can be expected. Nucleotides with phosphate groups show better peak shapes when compared to the regular column hardware. The **metal-free** column hardware is very suitable for highly sensitive analyses using LC/MS.


## IP-RP – Expert Tips: System and Column hardware

## Influence of system and column hardware on the analysis of nucleotides




1 ATP



2 ADP





3 AMP

Ordinary HPLC system  
with standard columnOrdinary HPLC system  
with metal-free column

Column: YMC-Triart C18 (3  $\mu$ m, 12 nm) 50 x 2.1 mm ID  
 Part Nos: TA12S03-05Q1PT (standard hardware)  
 TA12S03-05Q1PTP (metal-free hardware)  
 Eluent: 5 mM HCOONH<sub>4</sub>  
 Flow rate: 0.21 mL/min  
 Temperature: 25°C  
 Detection: UV at 265 nm  
 Injection: 1  $\mu$ L (10  $\mu$ g/mL)


"Non-metal" HPLC system: PEEK sample loop, PEEK injector port, and PEEK tubing are used.






ATP peak is detected, and peak shape of ADP is improved as a result of using the metal-free column.



Ordinary HPLC system  
with metal-free column



**"Non-metal" HPLC system\***  
with metal-free column




Peak shape is greatly improved as a result of using "non-metal" HPLC system

## IP-RP – Expert Tips: Ion-pairing salts

## Effect of composition and salt concentration of ion-pairing mobile phase on the separation and signal intensity

## Comparison of separation and ESI-MS signal intensity using different ion-pairing buffers and organic solvents



The mobile phase composition has different effects on the separation and signal intensity in electrospray ionisation mass spectrometry (ESI-MS) of oligonucleotides. Using different gradient conditions, acceptable retention and resolution can be achieved (upper UV chromatograms; red trace) for each separation by optimising the gradient slope of the organic solvent regardless of the type of mobile phase. The ESI-MS intensity is significantly influenced by

the type and concentration of ion-pairing buffer as shown in the lower MS chromatograms (orange trace). HFIP-TEA buffer/methanol systems provide the maximum MS intensity. Enhanced retention and MS intensity are obtained using 10 mM DBAA buffer compared to 10 mM TEAA buffer, and the lower DBAA concentration results in approximately 1.5–3 times increase in the intensity without any change in the concentration of organic solvent.




AEX



## AEX – BioPro IEX Series

## Features


- Porous or non-porous hydrophilic polymers
- High recovery of oligonucleotides
- Very high resolution
- Low nonspecific adsorption
- Excellent reproducibility



**strong anion  
exchanger**


|                               | BioPro IEX QA                           |
|-------------------------------|-----------------------------------------|
| Matrix                        | hydrophilic polymer (polymethacrylate)  |
| Particle size / $\mu\text{m}$ | 5                                       |
| Pore size / nm                | 100                                     |
| Charged group                 | $-\text{CH}_2\text{N}^+(\text{CH}_3)_3$ |
| Counter ion                   | $\text{Cl}^-$                           |
| Available pH range            | 2.0–12.0                                |
| Temperature range             | 4–60°C                                  |
| Pressure limit                | 2.5–3.5 MPa (360–510 psi)               |
| Column hardware               | PEEK                                    |

Also available in 10, 20, 30 or 75  $\mu\text{m}$  for preparative scale



Porous polymer beads

|                               | BioPro IEX QF                                                                      |
|-------------------------------|------------------------------------------------------------------------------------|
| Matrix                        | hydrophilic polymer (polymethacrylate)                                             |
| Particle size / $\mu\text{m}$ | 3, 5                                                                               |
| Pore size / nm                | non-porous                                                                         |
| Charged group                 | $-\text{CH}_2\text{N}^+(\text{CH}_3)_3$                                            |
| Counter ion                   | $\text{Cl}^-$                                                                      |
| Available pH range            | 2.0–12.0                                                                           |
| Temperature range             | 4–60°C                                                                             |
| Pressure limit                | 3 $\mu\text{m}$ : 25 MPa (3,625 psi)<br>5 $\mu\text{m}$ : 6–12 MPa (870–1,740 psi) |
| Column hardware               | PEEK                                                                               |




Non-porous polymer beads

YMC's anion exchanger (AEX) columns of the BioPro IEX series are available with strong exchanger modification, based on 5  $\mu\text{m}$  porous (QA columns) and on 3 or 5  $\mu\text{m}$  non-porous (QF columns) hydrophilic polymer beads.

The porous materials offer excellent binding capacity with exceptionally high efficiency and low operating pressure, whilst the non-porous particles offer high efficiency, very high resolution and low operating pressures.

## Column efficiency and loadability

### When to use porous and non-porous BioPro IEX




Column: BioPro IEX SF/SP  
 Eluent: A) 20 mM  $\text{NaH}_2\text{PO}_4$ - $\text{Na}_2\text{HPO}_4$  (pH 6.8)  
 B) 20 mM  $\text{NaH}_2\text{PO}_4$ - $\text{Na}_2\text{HPO}_4$  (pH 6.8) containing 0.5 M NaCl  
 Gradient: 0–100%B (0–60 min)  
 Flow rate: 0.5 mL/min  
 Temperature: 25°C

Detection: UV at 280 nm  
 Injection: 100 µL  
 Sample: 1. Ribonuclease A  
 2. Cytochrome c  
 3. Lysozyme

Non-porous BioPro IEX offers outstanding column efficiency at small amount of sample loading. Non-porous type of BioPro IEX columns are especially suitable for microscale analysis which requires higher resolution. Porous BioPro IEX maintains the good peak shape even when the loading amount increases. Porous type BioPro IEX columns with high capacity are useful for high-load analytical separations and laboratory-scale purification.

## AEX – BioPro IEX: Reproducibility

## Excellent batch-to-batch reproducibility



BioPro IEX columns exhibit excellent batch-to-batch reproducibility. All gel batches are inspected by rigorous quality control tests, and must meet the required criteria before release. BioPro IEX columns are the best choice for the quality control of biopharmaceuticals such as oligonucleotides or mAbs as in this example.

## Optimisation of oligonucleotide separations on ion exchange chromatography

Non-porous anion exchange column is generally suitable for analysis of oligonucleotides. For optimisation of single-stranded DNA and RNA of about 20 mer some conditions such as type of mobile phase and column temperature can be changed.


### 1 Improvement of peak tailing

#### Sample Group 1 (Phosphodiester oligonucleotides; PO)

|   |                             |                                                                                                                |
|---|-----------------------------|----------------------------------------------------------------------------------------------------------------|
| 1 | Single-stranded DNA (ssDNA) | 5'-TCATCACACTGAATACCAAT-3' (DNA 20 mer)                                                                        |
| 2 |                             | 5'-GTCATCACACTGAATACCAAT-3' (DNA 21 mer)                                                                       |
| 3 | Single-stranded RNA (ssRNA) | 5'-UCAUCACACUGAAUACCAU-3' (RNA 20 mer)                                                                         |
| 4 |                             | 5'-GUCAUCACACUGAAUACCAU-3' (RNA 21 mer)                                                                        |
| 5 |                             | 5'-U(M)C(M)A(M)U(M)C(M)A(M)C(M)A(M)C(M)U(M)G(M)A(M)A(M)U(M)A(M)C(M)C(M)A(M)A(M)U(M)-3' (2'-OMe RNA 20 mer)     |
| 6 |                             | 5'-G(M)U(M)C(M)A(M)U(M)C(M)A(M)C(M)A(M)C(M)U(M)G(M)A(M)A(M)U(M)A(M)C(M)C(M)A(M)A(M)U(M)-3' (2'-OMe RNA 21 mer) |

N(M)=2'-OMe RNA

By changing the buffer from 20 mM Tris-HCl (pH 8.1) to 10 mM NaOH, the tailing factor for an oligonucleotide is reduced. Furthermore, the peak tailing is further suppressed when NaClO<sub>4</sub> was added to 10 mM NaOH instead of NaCl.



### 2 Improvement of carryover

When the initial gradient concentration of NaCl is low (ex. 50 mM), carryover is observed. By increasing the initial gradient concentration of NaCl up to 400 mM, carryover can be avoided with good reproducibility.




## AEX Expert Tips

### 3 Improvement of ssDNA separation with single-base differences (differing in the type of base of 5'end of DNA 21mer)

When ssDNAs with single-base differences (differing in the type of base of 5'end of DNA 21mer) is analysed under the described condition, all of peak separations got worse. By using a shallower gradient, better separations could be achieved.

|   |                     |                                          |
|---|---------------------|------------------------------------------|
| 1 | Single-stranded DNA | 5'-TCATCACACTGAATACCAAT-3' (DNA 20 MER)  |
| 2 |                     | 5'-GTCATCACACTGAATACCAAT-3' (DNA 21 MER) |
| 7 |                     | 5'-ATCATCACACTGAATACCAAT-3' (DNA 21 MER) |
| 8 |                     | 5'-CTCATCACACTGAATACCAAT-3' (DNA 21 MER) |
| 9 |                     | 5'-TTCATCACACTGAATACCAAT-3' (DNA 21 MER) |



\*base of 5'end of DNA 21 mer



### 4 Improvement of the separation of phosphorothioate oligonucleotides with single-base differences in length

Since acidity of all PS is much higher than all PO, a higher salt concentration is required for elution. The peak of all PS is much broader because it is thought that all PS contains as many as 219 (524,288) stereoisomers. A steeper gradient curve, increasing column temperature and adding organic solvent can improve peak separation. However, increase of organic solvent ratio gave little improvement in peak separation.


### Sample Group 2 (Phosphorothioate oligonucleotides; PS)

|    |                     |                                                                     |
|----|---------------------|---------------------------------------------------------------------|
| 10 | Single-stranded DNA | 5'-T~T~T~T~T~T~T~T~T~T~T~T-3' (DNA 15 mer 12PS, 2PO)                |
| 11 |                     | 5'-T~T~T~T~T~T~T~T~T~T~T~T-3' (DNA 15 mer 13PS, 1PO)                |
| 12 |                     | 5'-T~T~T~T~T~T~T~T~T~T~T~T-3' (DNA 15 mer All PS)                   |
| 13 | Single-stranded RNA | 5'-U~C~A~U~C~A~C~A~C~U~G~A~A~U~A~C~C~A~A~U-3' (RNA 20 mer All PS)   |
| 14 |                     | 5'-G~U~C~A~U~C~A~C~A~C~U~G~A~A~U~A~C~C~A~A~U-3' (RNA 21 mer All PS) |

~ = Phosphorothioated







## Ordering information



## IP-RP – Ordering information

## 1.9 µm UHPLC columns

| Phase              | Column ID [mm] | Column length [mm] |                |                |                |                | Guard cartridges* with 5 mm length |
|--------------------|----------------|--------------------|----------------|----------------|----------------|----------------|------------------------------------|
|                    |                | 30                 | 50             | 75             | 100            | 150            | (pack of 3)                        |
| YMC-Triart C18     | 2.0            | TA12SP9-0302PT     | TA12SP9-0502PT | TA12SP9-L502PT | TA12SP9-1002PT | TA12SP9-1502PT | TA12SP9-E5Q1CC**                   |
|                    | 2.1            | TA12SP9-03Q1PT     | TA12SP9-05Q1PT | TA12SP9-L5Q1PT | TA12SP9-10Q1PT | TA12SP9-15Q1PT | TA12SP9-E5Q1CC**                   |
|                    | 3.0            | —                  | TA12SP9-0503PT | TA12SP9-L503PT | TA12SP9-1003PT | TA12SP9-1503PT | TA12SP9-E503CC                     |
| YMC-Triart Bio C18 | 2.0            | TA30SP9-0302PT     | TA30SP9-0502PT | TA30SP9-L502PT | TA30SP9-1002PT | TA30SP9-1502PT | TA30SP9-E5Q1CC**                   |
|                    | 2.1            | TA30SP9-03Q1PT     | TA30SP9-05Q1PT | TA30SP9-L5Q1PT | TA30SP9-10Q1PT | TA30SP9-15Q1PT | TA30SP9-E5Q1CC**                   |
|                    | 3.0            | —                  | TA30SP9-0503PT | TA30SP9-L503PT | TA30SP9-1003PT | TA30SP9-1503PT | TA30SP9-E503CC                     |
| YMC-Triart C8      | 2.0            | T012SP9-0302PT     | T012SP9-0502PT | T012SP9-L502PT | T012SP9-1002PT | T012SP9-1502PT | T012SP9-E5Q1CC**                   |
|                    | 2.1            | T012SP9-03Q1PT     | T012SP9-05Q1PT | T012SP9-L5Q1PT | T012SP9-10Q1PT | T012SP9-15Q1PT | T012SP9-E5Q1CC**                   |
|                    | 3.0            | —                  | T012SP9-0503PT | T012SP9-L503PT | T012SP9-1003PT | T012SP9-1503PT | T012SP9-E503CC                     |
| YMC-Triart Bio C4  | 2.0            | TB30SP9-0302PT     | TB30SP9-0502PT | TB30SP9-L502PT | TB30SP9-1002PT | TB30SP9-1502PT | TB30SP9-E5Q1CC**                   |
|                    | 2.1            | TB30SP9-03Q1PT     | TB30SP9-05Q1PT | TB30SP9-L5Q1PT | TB30SP9-10Q1PT | TB30SP9-15Q1PT | TB30SP9-E5Q1CC**                   |
|                    | 3.0            | —                  | TB30SP9-0503PT | TB30SP9-L503PT | TB30SP9-1003PT | TB30SP9-1503PT | TB30SP9-E503CC                     |

\*Guard cartridge holder required, part no. XPCUHHP

\*\*Guard cartridge: 2.1 mm ID

## 1.9 µm metal-free UHPLC columns

| Phase              | Column ID [mm] | Column length [mm] |                 |                 |
|--------------------|----------------|--------------------|-----------------|-----------------|
|                    |                | 50                 | 100             | 150             |
| YMC-Triart C18     | 2.1            | TA12SP9-05Q1PTP    | TA12SP9-10Q1PTP | TA12SP9-15Q1PTP |
| YMC-Triart Bio C18 | 2.1            | TA30SP9-05Q1PTP    | TA30SP9-10Q1PTP | TA30SP9-15Q1PTP |
| YMC-Triart C8      | 2.1            | T012SP9-05Q1PTP    | T012SP9-10Q1PTP | T012SP9-15Q1PTP |
| YMC-Triart Bio C4  | 2.1            | TB30SP9-05Q1PTP    | TB30SP9-10Q1PTP | TB30SP9-15Q1PTP |

Special column connectors required.

## Column connectors for metal-free (U)HPLC columns

| Recommendation  | ✓ ✓                                                                                 |                                                                                     | ✓                                                                                    |                                                                                       |
|-----------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| Ferrule         | no                                                                                  |                                                                                     | replaceable                                                                          |                                                                                       |
| Product         | MarvelX™                                                                            | MarvelXACT™                                                                         | Handy connector 2                                                                    | Hand-tight EXP® fitting                                                               |
| Manufacturer    | IDEX Health & Science LLC                                                           | IDEX Health & Science LLC                                                           | YMC Co., Ltd.                                                                        | Optimize Technologies, Inc.                                                           |
|                 |  |  |  |  |
| Pressure rating | 131 MPa / 1,310 bar                                                                 | 131 MPa / 1,310 bar                                                                 | 42 MPa / 420 bar                                                                     | 137 MPa / 1,370 bar                                                                   |
| Product code    | e.g. UPFP-6050250                                                                   | e.g. UPFP-YM7050250                                                                 | XRP0204                                                                              | XRHTF-01                                                                              |

MarvelX (ACT) is a registered trademark of IDEX Health &amp; Science LLC - EXP® is a registered trademark of Optimize Technologies, Inc.

## 3 µm HPLC columns

| Phase              | Column ID [mm] | Column length [mm] |                 |                 |                 |                 |                 | Guard cartridges* with 10 mm length |
|--------------------|----------------|--------------------|-----------------|-----------------|-----------------|-----------------|-----------------|-------------------------------------|
|                    |                | 30/33              | 50              | 75              | 100             | 150             | 250             |                                     |
|                    |                |                    |                 |                 |                 |                 |                 | (pack of 5)                         |
| YMC-Triart C18     | 2.1            | TA12S03-H3Q1PTH    | TA12S03-05Q1PTH | TA12S03-L5Q1PTH | TA12S03-10Q1PTH | TA12S03-15Q1PTH | –               | TA12S03-01Q1GC                      |
|                    | 3.0            | –                  | TA12S03-05Q3PTH | TA12S03-L5Q3PTH | TA12S03-10Q3PTH | TA12S03-15Q3PTH | –               | TA12S03-01Q3GC                      |
|                    | 4.6            | TA12S03-H346PTH    | TA12S03-0546PTH | TA12S03-L546PTH | TA12S03-1046PTH | TA12S03-1546PTH | TA12S03-2546PTH | TA12S03-01Q4GC                      |
| YMC-Triart Bio C18 | 2.1            | TA30S03-H3Q1PTH    | TA30S03-05Q1PTH | TA30S03-L5Q1PTH | TA30S03-10Q1PTH | TA30S03-15Q1PTH | –               | TA30S03-01Q1GC                      |
|                    | 3.0            | –                  | TA30S03-05Q3PTH | TA30S03-L5Q3PTH | TA30S03-10Q3PTH | TA30S03-15Q3PTH | –               | TA30S03-01Q3GC                      |
|                    | 4.6            | TA30S03-H346PTH    | TA30S03-0546PTH | TA30S03-L546PTH | TA30S03-1046PTH | TA30S03-1546PTH | TA30S03-2546PTH | TA30S03-01Q4GC                      |
| YMC-Triart C8      | 2.1            | T012S03-H3Q1PTH    | T012S03-05Q1PTH | T012S03-L5Q1PTH | T012S03-10Q1PTH | T012S03-15Q1PTH | –               | T012S03-01Q1GC                      |
|                    | 3.0            | –                  | T012S03-05Q3PTH | T012S03-L5Q3PTH | T012S03-10Q3PTH | T012S03-15Q3PTH | –               | T012S03-01Q3GC                      |
|                    | 4.6            | T012S03-H346PTH    | T012S03-0546PTH | T012S03-L546PTH | T012S03-1046PTH | T012S03-1546PTH | T012S03-2546PTH | T012S03-01Q4GC                      |
| YMC-Triart Bio C4  | 2.1            | TB30S03-03Q1PTH    | TB30S03-05Q1PTH | TB30S03-L5Q1PTH | TB30S03-10Q1PTH | TB30S03-15Q1PTH | –               | TB30S03-01Q1GC                      |
|                    | 3.0            | –                  | TB30S03-05Q3PTH | TB30S03-L5Q3PTH | TB30S03-10Q3PTH | TB30S03-15Q3PTH | –               | TB30S03-01Q3GC                      |
|                    | 4.6            | TB30S03-0346PTH    | TB30S03-0546PTH | TB30S03-L546PTH | TB30S03-1046PTH | TB30S03-1546PTH | TB30S03-2546PTH | TB30S03-01Q4GC                      |
| Hydrosphere C18    | 2.1            | HS12S03-03Q1WT     | HS12S03-05Q1WT  | HS12S03-L5Q1WT  | HS12S03-10Q1WT  | HS12S03-15Q1WT  | HS12S03-25Q1WT  | HS12S03-01Q1GC                      |
|                    | 3.0            | HS12S03-03Q3WT     | HS12S03-05Q3WT  | HS12S03-L5Q3WT  | HS12S03-10Q3WT  | HS12S03-15Q3WT  | HS12S03-25Q3WT  | HS12S03-01Q3GC                      |
|                    | 4.6            | HS12S03-0346WT     | HS12S03-0546WT  | HS12S03-L546WT  | HS12S03-1046WT  | HS12S03-1546WT  | HS12S03-2546WT  | HS12S03-01Q4GC                      |

\*Guard cartridge holder required, part no. XPGCH-Q1

## 3 µm metal-free HPLC columns

| Phase              | Column ID [mm] | Column length [mm] |                 |                 |
|--------------------|----------------|--------------------|-----------------|-----------------|
|                    |                | 50                 | 100             | 150             |
| YMC-Triart C18     | 2.1            | TA12S03-05Q1PTP    | TA12S03-10Q1PTP | TA12S03-15Q1PTP |
|                    | 4.6            | TA12S03-0546PTP    | TA12S03-1046PTP | TA12S03-1546PTP |
| YMC-Triart Bio C18 | 2.1            | TA30S03-05Q1PTP    | TA30S03-10Q1PTP | TA30S03-15Q1PTP |
|                    | 4.6            | TA30S03-0546PTP    | TA30S03-1046PTP | TA30S03-1546PTP |
| YMC-Triart C8      | 2.1            | T012S03-05Q1PTP    | T012S03-10Q1PTP | T012S03-15Q1PTP |
|                    | 4.6            | T012S03-0546PTP    | T012S03-1046PTP | T012S03-1546PTP |
| YMC-Triart Bio C4  | 2.1            | TB30S03-05Q1PTP    | TB30S03-10Q1PTP | TB30S03-15Q1PTP |
|                    | 4.6            | TB30S03-0546PTP    | TB30S03-1046PTP | TB30S03-1546PTP |

Special column connectors required.

## IP-RP – Ordering information

## 5 µm HPLC columns

| Phase              | Column ID [mm] | Column length [mm] |                 |                 |                 |                 |                 | Guard cartridges* with 10 mm length |
|--------------------|----------------|--------------------|-----------------|-----------------|-----------------|-----------------|-----------------|-------------------------------------|
|                    |                | 30/33              | 50              | 75              | 100             | 150             | 250             |                                     |
| YMC-Triart C18     | 2.1            | TA12S05-H3Q1PTH    | TA12S05-05Q1PTH | TA12S05-L5Q1PTH | TA12S05-10Q1PTH | TA12S05-15Q1PTH | –               | TA12S05-01Q1GC                      |
|                    | 3.0            | –                  | TA12S05-05Q3PTH | TA12S05-L5Q3PTH | TA12S05-10Q3PTH | TA12S05-15Q3PTH | –               | TA12S05-01Q3GC                      |
|                    | 4.6            | TA12S05-H346PTH    | TA12S05-0546PTH | TA12S05-L546PTH | TA12S05-1046PTH | TA12S05-1546PTH | TA12S05-2546PTH | TA12S05-01Q4GC                      |
| YMC-Triart Bio C18 | 2.1            | TA30S05-H3Q1PTH    | TA30S05-05Q1PTH | TA30S05-L5Q1PTH | TA30S05-10Q1PTH | TA30S05-15Q1PTH | –               | TA30S05-01Q1GC                      |
|                    | 3.0            | –                  | TA30S05-05Q3PTH | TA30S05-L5Q3PTH | TA30S05-10Q3PTH | TA30S05-15Q3PTH | –               | TA30S05-01Q3GC                      |
|                    | 4.6            | TA30S05-H346PTH    | TA30S05-0546PTH | TA30S05-L546PTH | TA30S05-1046PTH | TA30S05-1546PTH | TA30S05-2546PTH | TA30S05-01Q4GC                      |
| YMC-Triart C8      | 2.1            | T012S05-H3Q1PTH    | T012S05-05Q1PTH | T012S05-L5Q1PTH | T012S05-10Q1PTH | T012S05-15Q1PTH | –               | T012S05-01Q1GC                      |
|                    | 3.0            | –                  | T012S05-05Q3PTH | T012S05-L5Q3PTH | T012S05-10Q3PTH | T012S05-15Q3PTH | –               | T012S05-01Q3GC                      |
|                    | 4.6            | T012S05-H346PTH    | T012S05-0546PTH | T012S05-L546PTH | T012S05-1046PTH | T012S05-1546PTH | T012S05-2546PTH | T012S05-01Q4GC                      |
| YMC-Triart Bio C4  | 2.1            | TB30S05-H3Q1PTH    | TB30S05-05Q1PTH | TB30S05-L5Q1PTH | TB30S05-10Q1PTH | TB30S05-15Q1PTH | –               | TB30S05-01Q1GC                      |
|                    | 3.0            | –                  | TB30S05-05Q3PTH | TB30S05-L5Q3PTH | TB30S05-10Q3PTH | TB30S05-15Q3PTH | –               | TB30S05-01Q3GC                      |
|                    | 4.6            | TB30S05-H346PTH    | TB30S05-0546PTH | TB30S05-L546PTH | TB30S05-1046PTH | TB30S05-1546PTH | TB30S05-2546PTH | TB30S05-01Q4GC                      |
| Hydrosphere C18    | 2.1            | HS12S05-03Q1WT     | HS12S05-05Q1WT  | HS12S05-L5Q1WT  | HS12S05-10Q1WT  | HS12S05-15Q1WT  | HS12S05-25Q1WT  | HS12S05-01Q1GC                      |
|                    | 3.0            | HS12S05-03Q3WT     | HS12S05-05Q3WT  | HS12S05-L5Q3WT  | HS12S05-10Q3WT  | HS12S05-15Q3WT  | HS12S05-25Q3WT  | HS12S05-01Q3GC                      |
|                    | 4.6            | HS12S05-0346WT     | HS12S05-0546WT  | HS12S05-L546WT  | HS12S05-1046WT  | HS12S05-1546WT  | HS12S05-2546WT  | HS12S05-01Q4GC                      |

\*Guard cartridge holder required, part no. XPGCH-Q1

## 5 µm metal-free HPLC columns

| Phase              | Column ID [mm] | Column length [mm] |                 |                 |
|--------------------|----------------|--------------------|-----------------|-----------------|
|                    |                | 50                 | 100             | 150             |
| YMC-Triart C18     | 2.1            | TA12S05-05Q1PTP    | TA12S05-10Q1PTP | TA12S05-15Q1PTP |
|                    | 4.6            | TA12S05-0546PTP    | TA12S05-1046PTP | TA12S05-1546PTP |
| YMC-Triart Bio C18 | 2.1            | TA30S05-05Q1PTP    | TA30S05-10Q1PTP | TA30S05-15Q1PTP |
|                    | 4.6            | TA30S05-0546PTP    | TA30S05-1046PTP | TA30S05-1546PTP |
| YMC-Triart C8      | 2.1            | T012S05-05Q1PTP    | T012S05-10Q1PTP | T012S05-15Q1PTP |
|                    | 4.6            | T012S05-0546PTP    | T012S05-1046PTP | T012S05-1546PTP |
| YMC-Triart Bio C4  | 2.1            | TB30S05-05Q1PTP    | TB30S05-10Q1PTP | TB30S05-15Q1PTP |
|                    | 4.6            | TB30S05-0546PTP    | TB30S05-1046PTP | TB30S05-1546PTP |

Special column connectors required.

## 5 µm YMC-Actus high-throughput semipreparative columns

| Phase              | Column ID [mm] | Column length [mm] |                |                |                |                | Guard cartridges* with 10 mm length |
|--------------------|----------------|--------------------|----------------|----------------|----------------|----------------|-------------------------------------|
|                    |                | 50                 | 75             | 100            | 150            | 250            |                                     |
| YMC-Triart C18     | 20             | TA12S05-0520WX     | TA12S05-L520WX | TA12S05-1020WX | TA12S05-1520WX | TA12S05-2520WX | TA12S05-0120CCN                     |
|                    | 30             | TA12S05-0530WX     | TA12S05-L530WX | TA12S05-1030WX | TA12S05-1530WX | TA12S05-2530WX | TA12S05-0130CCN                     |
|                    | 50             | TA12S05-0553DX     | –              | TA12S05-1053DX | TA12S05-1553DX | TA12S05-2553DX | TA12S05-0553DXG**                   |
| YMC-Triart Bio C18 | 20             | TA30S05-0520WX     | TA30S05-L520WX | TA30S05-1020WX | TA30S05-1520WX | TA30S05-2520WX | TA30S05-0120CCN                     |
|                    | 30             | TA30S05-0530WX     | TA30S05-L530WX | TA30S05-1030WX | TA30S05-1530WX | TA30S05-2530WX | TA30S05-0130CCN                     |
|                    | 50             | TA30S05-0553DX     | –              | TA30S05-1053DX | TA30S05-1553DX | TA30S05-2553DX | TA30S05-0553DXG**                   |
| YMC-Triart C8      | 20             | T012S05-0520WX     | T012S05-L520WX | T012S05-1020WX | T012S05-1520WX | T012S05-2520WX | T012S05-0120CCN                     |
|                    | 30             | T012S05-0530WX     | T012S05-L530WX | T012S05-1030WX | T012S05-1530WX | T012S05-2530WX | T012S05-0130CCN                     |
|                    | 50             | T012S05-0553DX     | –              | T012S05-1053DX | T012S05-1553DX | T012S05-2553DX | T012S05-0553DXG**                   |
| YMC-Triart Bio C4  | 20             | TB30S05-0520WX     | TB30S05-L520WX | TB30S05-1020WX | TB30S05-1520WX | TB30S05-2520WX | TB30S05-0120CCN                     |
|                    | 30             | TB30S05-0530WX     | TB30S05-L530WX | TB30S05-1030WX | TB30S05-1530WX | TB30S05-2530WX | TB30S05-0130CCN                     |
|                    | 50             | TB30S05-0553DX     | –              | TB30S05-1053DX | TB30S05-1553DX | TB30S05-2553DX | TB30S05-0553DXG**                   |
| Hydrosphere C18    | 20             | HS12S05-0520WX     | HS12S05-L520WX | HS12S05-1020WX | HS12S05-1520WX | HS12S05-2520WX | HS12S05-0120CCN                     |
|                    | 30             | HS12S05-0530WX     | HS12S05-L530WX | HS12S05-1030WX | HS12S05-1530WX | HS12S05-2530WX | HS12S05-0130CCN                     |
|                    | 50             | HS12S05-0553DX     | –              | HS12S05-1053DX | HS12S05-1553DX | HS12S05-2553DX | –                                   |

\*Guard cartridge holder required, part no. XPGHF2P20ID (20 mm ID)  
 XPGHF2P30ID (30 mm ID)  
 no holder required for 50 mm



## AEX – Ordering information

**3 µm non-porous analytical columns, PEEK hardware (max. pressure 250 bar)**

| Phase         | Column ID [mm] | Column length [mm] |                |                |                | Precolumn filter 2 µm* |
|---------------|----------------|--------------------|----------------|----------------|----------------|------------------------|
|               |                | 30 (250 bar)       | 50 (250 bar)   | 100 (250 bar)  | 150 (250 bar)  | (pack of 5)            |
| BioPro IEX QF | 4.6            | QF00S03-0346WP     | QF00S03-0546WP | QF00S03-1046WP | QF00S03-1546WP | XRPRCP25               |

**5 µm non-porous analytical columns, PEEK hardware (max. pressure 60–120 bar)**

| Phase         | Column ID [mm] | Column length [mm] |                |                |                | Precolumn filter 2 µm* |
|---------------|----------------|--------------------|----------------|----------------|----------------|------------------------|
|               |                | 30 (60 bar)        | 50 (100 bar)   | 100 (120 bar)  | 150 (120 bar)  | (pack of 5)            |
| BioPro IEX QF | 4.6            | QF00S05-0346WP     | QF00S05-0546WP | QF00S05-1046WP | QF00S05-1546WP | XRPRCP25               |

**5 µm porous analytical columns, PEEK hardware (max. pressure 25–35 bar)**

| Phase         | Column ID [mm] | Column length [mm] |                |                | Precolumn filter 2 µm* |
|---------------|----------------|--------------------|----------------|----------------|------------------------|
|               |                | 30 (25 bar)        | 50 (30 bar)    | 100 (35 bar)   | (pack of 5)            |
| BioPro IEX QA | 4.6            | QAA0S05-0346WP     | QAA0S05-0546WP | QAA0S05-1046WP | XRPRCP25               |

\* Holder required, part no. XRPRCP02

**6 µm non-porous semiprep. columns, stainless steel hardware (max. pressure 30–90 bar)**

| Phase         | Column ID [mm] | Column length [mm]                                 |  |
|---------------|----------------|----------------------------------------------------|--|
|               |                | 100                                                |  |
| BioPro IEX QF | 10<br>20<br>30 | QF00S06-1010WT<br>QF00S06-1020WT<br>QF00S06-1030WT |  |

**6 µm porous semiprep. columns, stainless steel hardware (max. pressure 40 bar)**

| Phase         | Column ID [mm] | Column length [mm]               |  |
|---------------|----------------|----------------------------------|--|
|               |                | 100                              |  |
| BioPro IEX QA | 10<br>20       | QAA0S06-1010WT<br>QAA0S06-1020WT |  |

Other dimensions on demand



Your local distributor:

**YMC** America, Inc.  
8 Charlestown Street  
Devens, MA 01434  
Phone +1 978 487 1100  
Email: [info@ymcamerica.com](mailto:info@ymcamerica.com)  
[www.ymcamerica.com](http://www.ymcamerica.com)

**YMC** Europe GmbH  
Schöttmannshof 19  
D-46539 Dinslaken  
Germany  
Phone +49(0)2064 427-0, Fax +49(0)2064 427-222  
[www.ymc.eu](http://www.ymc.eu)

**YMC** Schweiz GmbH  
Im Wasenboden 8  
4056 Basel  
Switzerland  
Phone +41(0)61 561 8050, Fax +41(0)61 561 8059  
[www.ymc-schweiz.ch](http://www.ymc-schweiz.ch)

**YMC** CO., LTD.  
YMC Karasuma-Gojo Bld. 284 Daigo-cho,  
Karasuma Nishiiru Gojo-dori Shimogyo-ku,  
Kyoto 600-8106 Japan  
Phone +81(0)75 342 4515, Fax +81(0)75 342 4550  
[www.ymc.co.jp](http://www.ymc.co.jp)