


# YMC Glass Column Packing Service

## Save time with pre-packed YMC glass columns!

Packing columns can be a time-consuming task as well as a drain on resources.

Use YMC's packing expertise to your advantage and allow us to pack your glass columns for you!



## Available Services:

*Packing of glass columns according to your IEX, Affinity, SEC, HIC, or RP applications*

- Inner diameter of 5 - 80mm
- Length of up to 1,000mm
- Customized dimensions possible

*We can pack new columns or refill your existing columns*

- Order a new YMC glass column with packing services, and we will ship it to you packed as directed and ready for use!
- Send us your existing YMC glass column to be re-filled. Additional maintenance services, such as changing frits & o-rings, can be performed at the same time.

*Choose stationary phases from YMC or other suppliers*

Of course we recommend YMC packing materials, but we understand that sometimes you have to follow existing specifications. In these cases, we are still here to help! Here are the options:

- Order new YMC resins to be packed in your glass column
- Send us the required quantity from your existing YMC inventory
- Send us the resin to be packed from an alternate supplier



Complete the form on reverse to receive a quote.



# Glass Column Packing Service

For help in selecting your laboratory or pilot scale glass column, please complete this questionnaire and return it to your rep or to [info@ymcamerica.com](mailto:info@ymcamerica.com).

1. Inner diameter (mm): \_\_\_\_\_
2. Minimum bed length (mm): \_\_\_\_\_
3. Maximum bed length (mm): \_\_\_\_\_
4. Bed volume (mL): \_\_\_\_\_
5. Maximum pressure (bar): \_\_\_\_\_
6. Biocompatibility: Yes No
7. Packing material: \_\_\_\_\_
8. Particle size of packing material: \_\_\_\_\_
9. Solvents: \_\_\_\_\_
10. Heating/cooling jacket: Yes No
11. Quotation for packing adapter required: Yes No
12. Regulatory requirements: \_\_\_\_\_
13. Additional notes: \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_